



Test Update 711

Posted Date 04/01/2020
Effective Date 04/09/2020
Update Type [New Tests](#)
CPT Code 80299, 82397

Order Code: FGAGA

Reference Laboratory: Mayo Medical Laboratories

Effective April 9, 2020 MLabs will offer Golimumab and Anti-Golimumab Ab (FGAGA) testing via Mayo Medical Laboratories.

Specimen Collection: Collect approximately 6 mL of whole blood in a red or serum separator tube (SST). Centrifuge the specimen within 45 minutes of collection. Aliquot the serum into a plastic vial and freeze. Ship the specimen frozen within 7 days of collection.

Methodology: Electrochemiluminescence immunoassay (ECLIA)

Analytic Time: Approx. 2 days, test performed on Tuesdays

Reference Range: See Below

Golimumab:

Quantitation Limit: <0.5 ug/mL

Results of 0.5 ug/mL or higher indicate detection of Golimumab

In the presence of serum anti-golimumab antibodies, the golimumab drug level reflects the antibody-unbound (free) fraction of golimumab in serum

Anti-Golimumab Antibody

2800 Plymouth Rd, Bldg 35 • Ann Arbor, MI 48109-2800 • 800.862.7284 • 734.936.0755 Fax • mlabs.umich.edu

Quantitation Limit: <20 ng/mL

Results of 20 or higher indicate detection of anti-Golimumab antibodies.

[View PDF](#)