



## Test Update 711

**Posted Date** 04/01/2020  
**Effective Date** 04/09/2020  
**Update Type** [New Tests](#)  
**CPT Code** 80299, 82397

**Order Code:** FGAGA

**Reference Laboratory:** Mayo Medical Laboratories

Effective April 9, 2020 MLabs will offer Golimumab and Anti-Golimumab Ab (FGAGA) testing via Mayo Medical Laboratories.

**Specimen Collection:** Collect approximately 6 mL of whole blood in a red or serum separator tube (SST). Centrifuge the specimen within 45 minutes of collection. Aliquot the serum into a plastic vial and freeze. Ship the specimen frozen within 7 days of collection.

**Methodology:** Electrochemiluminescence immunoassay (ECLIA)

**Analytic Time:** Approx. 2 days, test performed on Tuesdays

**Reference Range:** See Below

### **Golimumab:**

**Quantitation Limit:** <0.5 ug/mL

Results of 0.5 ug/mL or higher indicate detection of Golimumab

In the presence of serum anti-golimumab antibodies, the golimumab drug level reflects the antibody-unbound (free) fraction of golimumab in serum

### **Anti-Golimumab Antibody**

2800 Plymouth Rd, Bldg 35 • Ann Arbor, MI 48109-2800 • 800.862.7284 • 734.936.0755 Fax • [mlabs.umich.edu](http://mlabs.umich.edu)

**Quantitation Limit:** <20 ng/mL

Results of 20 or higher indicate detection of anti-Golimumab antibodies.

[View PDF](#)