



## **Test Update 853**

**Posted Date** 08/16/2023

**Effective Date** 08/14/2023

**Test Name** [Fungitell, BAL](#)

**Update Type** [Specimen collection or handling and/or storage changed](#)

**CPT Code** 87449

### **SPECIMEN HANDLING CHANGE**

#### **Fungitell, BAL**

Order Code: VFNFL

Fee Code: AA238 (CPT 87449)

Reference Laboratory: Viracor 351700

Viracor has updated the preferred specimen storage and transport to **room (ambient) temperature** for the Fungitell assays effective August 14, 2023.

Collection Instructions BAL: Collect 1-3 mL in a sterile, screw top tube. Store and send at room temperature (preferred), refrigerated, or frozen. Do not ship samples in universal transport media.

[View PDF](#)