



## **Test Update 685**

**Posted Date** 10/02/2019  
**Effective Date** 10/24/2019  
**Update Type** [New Tests](#)  
**CPT Code** 80194

Order Code: QUIN

Referral Laboratory: Mayo Medical Laboratories

Collection Instructions: Collect approximately 1 mL of whole blood in a red top vacutainer, centrifuge and aliquot the serum into a plastic vial within 2 hours of collection.

Methodology: Kinetic Interaction of Microparticles in Solution (KIMS)

Reference Values: **Therapeutic:** 2.0-5.0 mcg/mL **Critical Value:**  $\geq 6.0$  mcg/mL

Set Up: Mon- Sat Analytic Time: 1 Day

NOTE:

### **TEST OBSOLETE**

#### **Quinidine**

Order Code: QUIND

CPT Code: 80194

Referral Laboratory: Mayo Medical Laboratories

Effective October 24, 2019, MLabs will no longer offer Quinidine (QUIND) testing due to Mayo Medical laboratories discontinuation. The test will be replaced by a new Quinidine with a change in methodology.

[View PDF](#)