



## **Test Update 734**

**Posted Date** 10/14/2020  
**Effective Date** 10/28/2020  
**Test Name** [Androstenediol Glucuronide](#)  
**Update Type** [Replacement Test](#)  
**CPT Code** 82154

### **TEST CHANGE**

#### **Androstenediol Glucuronide**

**Order Code:** MMLR  
**CPT Code:** 82154  
**Referral Laboratory:** Mayo Medical Laboratories via Esoterix

### **REPLACEMENT TEST**

#### **Androstenediol Glucuronide**

**Order Code:** MMLR  
**CPT Code:**  
**Referral Laboratory:** Mayo Medical Laboratories via Esoterix Endocrinology

Effective October 28, 2020 MLabs will discontinue the Androstenediol Glucuronide (MMLR) test per Mayo Medical Laboratories via Esoterix and will be replaced by Androstenediol Glucuronide (MMLR) test per Mayo Medical Laboratories via Esoterix Endocrinology

Collection Instruction: Please collect 5 mL of whole blood in a red top (preferred) or serum separator tube. Centrifuge and aliquot approximately 1 mL serum into a plastic vial and ship refrigerated to MLabs within 6 days of collection.

Methodology: High Pressure Liquid Chromatography and Tandem Mass Spectrometry (HPLC/MS-MS)

Analytic Time: Set up on Mondays, 8-15 days

Reference Values:

**Age Range****Range**

Prepubertal Children Not Established

N/A

Adult Males

112 - 1046  
ng/dL

Adult Females

11 - 249  
ng/dL

\*\*Occasionally, normal females with no evidence of hirsutism may have levels well beyond the normal range

[View PDF](#)