



## **Test Update 735**

**Posted Date** 10/28/2020  
**Effective Date** 10/28/2020  
**Update Type** [New Tests](#)  
**CPT Code** 83520

### **NEW TEST**

#### **Inhibin A & B, Tumor Marker**

**Order Code:** INHAB  
**CPT Code:** 83520  
**Referral Laboratory:** Mayo Medical Laboratories

Effective October 28, 2020, MLabs will offer Inhibin A & B, Tumor Marker (INHAB) via Mayo Medical Laboratories.

Collection Instructions: Collect approximately 2 mL of whole blood in a red top vacutainer or a serum separator tube (SST). Centrifuge and aliquot 1 mL of serum into a plastic vial. Ship the specimen refrigerated.

Methodology: ELISA

Analytic Time: 1 day

Reference Range: Age and sex specific ranges provided with report

Test Usage: As an aid in the diagnosis of granulosa cell tumors and mucinous epithelial ovarian tumors and for monitoring of patients with granulosa cell tumors and epithelial mucinous-type tumors of the ovary known to overexpress inhibin A and B

[View PDF](#)