



## **Test Update 735**

**Posted Date** 10/28/2020  
**Effective Date** 11/12/2020  
**Update Type** [New Tests](#)  
**CPT Code** 83652 x2

### **NEW TEST**

#### **DHR Flow, Blood**

**Order Code:** DHR  
**CPT Code:** 83652 x2  
**Referral Laboratory:** Mayo Medical Laboratories

Effective October 28, 2020, MLabs will offer DHR Flow, Blood (DHR) via Mayo Medical Laboratories.

Collection Instructions: Collect 5 mL of whole blood in a sodium heparin green top vacutainer. Ship the specimen ambient within 24 hours of collection.

Methodology: Flow Cytometry

Analytic Time: 3-4 days

Reference Range:

<b>Result Name</b>	<b>Unit</b>	<b>Cutoff for Defining Normal</b>
--------------------	-------------	-----------------------------------

%PMA ox-DHR+	%	≥ 95%
--------------	---	-------

MFI PMA ox-DHR+	MFI	≥ 60
-----------------	-----	------

% fMLP ox-DHR+	%	≥ 10%
----------------	---	-------

MFI fMLP ox-DHR+	MFI	≥ 2
------------------	-----	-----

2000 Plymouth Rd, LBN 735 • Ann Arbor, MI 48109-2800 • 800.862.7284 • 734.936.0755 Fax • [mlabs.umich.edu](http://mlabs.umich.edu)

Control % PMA ox-DHR+ %  $\geq 95\%$

Control MFI PMA ox-DHR+ MFI  $\geq 60$

Control % fMLP ox-DHR+ %  $\geq 10\%$

Control MFI fMLP ox-DHR+ MFI  $\geq 2$

[View PDF](#)