



## [Test Update 737](#)

**Posted Date** 11/11/2020  
**Effective Date** 11/12/2020  
**Update Type** [New Tests](#)  
**CPT Code** 86352

Effective November 12, 2020, DHR (DHR) testing will be available via Mayo Medical Laboratories.

Collection Instructions: Collect 5 mL of whole blood in a sodium heparin green top vacutainer from both a patient and unrelated healthy donor within an hour of the patient. Clearly label each vacutainer control or patient.

Methodology: Flow Cytometry

Analytic Time: M-F, 4 days turnaround time

[View PDF](#)