



[Test Update 738](#)

Posted Date 11/18/2020
Effective Date 12/07/2020
Update Type [New Tests](#)
CPT Code 86300

Breast Carcinoma-Associated Antigen (CA27.29)

Order Code: 2729
CPT Code: 86300

On December 7, 2020 MLabs will offer Breast Carcinoma-Associated Antigen (CA27.29) [2729] testing in-house.

Test Usage: The Breast Carcinoma-Associated Antigen (CA27.29) assay is intended for use as an aid in monitoring patients previously treated for Stage II or Stage III breast cancer. It is also intended for use as an aid in the management of breast cancer patients with metastatic disease by monitoring the progression or regression of disease in response to treatment.

Collection Instructions: Collect approximately 3 mL of whole blood in a serum separator tube (SST). Centrifuge and aliquot the serum into a plastic vial within 2 hours of collection. Ship the specimen frozen.

Methodology: Chemiluminescent Immunoassay

Analytic Time: 2-3 days

Reference Range: Age 18 and over 0-38 U/mL

[View PDF](#)