



## Test Update 739

**Posted Date** 12/02/2020

**Effective Date** 12/03/2020

**Update Type** [Replacement Test](#)

**CPT Code** 83520

### TEST OBSOLETE

#### Anti-Mullerian Hormone (AMH)

**Order Code:** AMHS

**CPT Code:** 83520

**Referral Laboratory:** Mayo Medical Laboratories

Effective December 3, 2020, Anti-Mullerian Hormone (AMHS) testing will be replaced by Anti-Mullerian Hormone (AMH1) per Mayo Medical Laboratories.

### REPLACEMENT TEST

#### Anti-Mullerian Hormone (AMH)

**Order Code:** AMH1

**CPT Code:** 83520

**Referral Laboratory:** Mayo Medical Laboratories

**Collection Instructions:** Collect approximately 3 mL of whole blood in a serum separator tube (SST) or red top vacutainer. Centrifuge and aliquot a minimum 1 mL serum into a plastic vial. Ship the specimen refrigerators (preferred), frozen or ambient.

**Methodology:** Electrochemiluminescence assay

**Analytic Time:** 1 day turnaround time

#### Reference Range:

2800 Plymouth Rd, Bldg 35 • Ann Arbor, MI 48109-2800 • 800.862.7284 • 734.936.0755 Fax • mlabs.umich.edu

**Males**      **Range ng/mL**

<2 years:      18-283 ng/mL

2-12 years:      8.9-109 ng/mL

>12 years:      <13 ng/mL

**Females**      **Range ng/mL**

<3 years:      0.11-4.2  
ng/mL

3-6 years:      0.21-4.9  
ng/mL

7-11 years:      0.36-5.9  
ng/mL

12-14 years:      0.49-6.9  
ng/mL

15-19 years:      0.62-7.8  
ng/mL

20-24 years:      1.2-12 ng/mL

25-29 years:      0.89-9.9  
ng/mL

30-34 years:      0.58-8.1  
ng/mL

35-39 years:      0.15-7.5  
ng/mL

40-44 years:      0.03-5.5  
ng/mL

45-50 years: <2.6 ng/mL

51-55 years: <0.88 ng/mL

>55 years: <0.03 ng/mL

[View PDF](#)