



## **Test Update 746**

**Posted Date** 01/20/2021  
**Effective Date** 01/20/2021  
**Update Type** [New Tests](#)  
**CPT Code** 83520

**Reference Laboratory:** Esoterix

Effective January 20, 2021, MLabs will offer Anti-cN-1A Ab (NT5c1A) [LCN1A] testing via Esoterix.

Collection Instructions: Collect approximately 5 mL of whole blood in a serum separator tube (SST). Centrifuge and aliquot serum into a plastic vial within one hour of collection. Ship the specimen refrigerated.

Methodology: Enzyme-linked immunosorbent assay (ELISA)

Test Analysis: 14-21 days

[View PDF](#)