



Test Update 762

Posted Date 05/12/2021
Effective Date 05/26/2021
Test Name [Plasminogen Activator Inhibitor Antigen, Plasma](#)
Update Type [Replacement Test](#)
CPT Code 85415

TEST DOWN

Plasminogen Activator Inhibitor Antigen, Plasma

Order Code: PAI1
CPT Code: 85415
Referral Laboratory: Mayo Medical Laboratory

TEMPORARY REPLACEMENT

PAI-1 Antigen

Order Code: FPAI1
CPT Code: 85415
Referral Laboratory: Mayo Medical Laboratory

Effective May 26, 2021, Plasminogen Activator Inhibitor Antigen, Plasma (PAI1) testing will be down and the recommended alternative is PAI-1 Antigen (FPAI1) via Mayo Medical Laboratories.

Collection Instructions: Collect approximately 2 mL of whole blood in a sodium citrate vacutainer. Centrifuge and aliquot the plasma in a plastic tube and ship the specimen frozen.

Methodology: Enzyme-Linked Immunosorbent Assay (ELISA)

Test Analysis: Set Up Friday, 4-8 days turnaround time

Reference Ranges: 4.0 - 43.0

