



## **Test Update 768**

**Posted Date** 07/14/2021

**Effective Date** 06/30/2021

**Test Name** [Connective Tissue Gene Tests \(CTGT\), saliva kit](#)

**Update Type** [New Tests](#)

### **NEW TEST**

#### **Connective Tissue Gene Tests (CTGT), saliva kit**

**Order Code:** CTSAL

**Referral Laboratory:** Connective Tissue Gene Tests

Effective June 30, 2021, MLabs will offer Connective Tissue Gene Tests (CTGT), saliva kit (CTSAL) via Connective Tissue Gene Tests.

Collection Instructions: Collect saliva in a saliva collection container. Ship the specimen at room temperature

Methodology: PCR

Analytic Time: 2-4 weeks

[View PDF](#)