



Test Update 775

Posted Date 09/02/2021

Effective Date 09/21/2021

Test Name [Inhibin A and B, Tumor Marker](#)

Update Type [Specimen collection or handling and/or storage changed](#)

CPT Code 83520

TEST CHANGE

Inhibin A and B, Tumor Marker

Order Code: INHAB

CPT Code: 83520

Referral Laboratory: Mayo Medical Laboratories

Effective September 21, 2021, MLabs has changed the specimen requirements and stability for Inhibin A and B, Tumor Marker (INHAB) per Mayo Medical Laboratories.

Current Specimen Required

Collection Tube/ Container:

Preferred Red Top

Acceptable Serum Gel

New Specimen Required

Collection Tube/ Container:

Preferred Serum Gel

Acceptable Red Top

**Current Specimen
Stability**

Specimen	Temperature	Time
Serum	Refrigerated	7 days
	Frozen	90 Days

**New Specimen
Stability**

Specimen	Temperature	Time
Serum	Refrigerated	14 Days
	Frozen	90 Days

[View PDF](#)