



## **Test Update 777**

**Posted Date** 09/30/2021  
**Effective Date** 09/27/2021  
**Test Name** [Lead Profile Occupational Exposure](#)  
**Update Type** [Replacement Test](#)  
**CPT Code** 83655, 84202

### **TEST CHANGE**

#### **Lead Profile Occupational Exposure**

**Order Code:** PBZP  
**CPT Code:** 83655, 84202  
**Referral Laboratory:** Mayo Medical Laboratories

### **REPLACEMENT TEST**

#### **Lead, Industrial Profile, Adults**

**Order Code:** ALEAD  
**CPT Code:** 83655; 84202  
**Referral Laboratory:** ARUP

Effective September 27, 2021, MLabs will temporarily send Lead Profile Occupational Exposure (PBZP) to ARUP while Mayo Medical Laboratories testing is down.

Collection Instructions: Collect 7 ml of whole blood in a Royal blue (K2EDTA or Na2EDTA) or tan (K2EDTA) vacutainer. Ship specimen refrigerated. Specimen must be collected from a patient that is 16 years or older.

Analytic Time: 1-5 days

Reference Range: Interpretive

[View PDF](#)