



Test Update 805

Posted Date 05/11/2022

Effective Date 05/10/2022

Test Name [Nicotine and Metabolites, Urine](#)

Update Type [Specimen collection or handling and/or storage changed](#)

SPECIMEN HANDLING CHANGE

Nicotine and Metabolites, Urine

Order Code: NICOU

Fee Code: 21609

Please note that effective May 10, 2022, the preferred specimen temperature for the Nicotine and Metabolites, Urine, assay has been updated from ambient to refrigerated.

Collection Instructions: Collect random urine specimen. Refrigerate up to 28 days or freeze for longer storage.

[View PDF](#)