



## **Test Update 816**

**Posted Date** 08/24/2022

**Effective Date** 09/30/2022

**Test Name** [Plasminogen Activator Inhibitor Antigen, Plasma](#)

**Update Type** [Test Resumed](#)

TEST RESUMED

Plasminogen Activator Inhibitor Antigen, Plasma

Order Code: PAI1

Fee Code: AA606

Reference Laboratory: Mayo PAI1

Mayo Clinic Laboratories Plasminogen Activator Inhibitor Antigen, Plasma testing will resume effective August 30, 2022. This test replaces Mayo test ID FPAI1 referred to Esoterix.

Collection Instructions: Collect specimen in light blue (sodium citrate) tube. Centrifuge for 10 minutes and carefully aliquot plasma into a plastic vial, being careful not to disturb cells. Re-centrifuge for 10 minutes and aliquot 1 mL (minimum 0.5 mL) of plasma into a plastic vial, leaving platelets at the bottom of the tube. Freeze.

[View PDF](#)