



## **Test Update 840**

**Posted Date** 04/19/2023

**Effective Date** 04/28/2023

**Test Name** [Entamoeba histolytica \(Amoeba\) Antibody](#)

**Update Type** [Specimen collection or handling and/or storage changed](#)

### **SPECIMEN COLLECTION AND HANDLING UPDATE**

#### **Entamoeba histolytica (Amoeba) Antibody**

Order Code: EHOLG

Fee Code: 32052 (CPT 86753)

Reference Laboratory: Mayo EHOLG

The specimen stability for Mayo Clinic Laboratories Entamoeba histolytica Antibody assay will be updated effective April 28, 2023.

Collection Instructions: Collect specimen in an SST or red top tube. Centrifuge, aliquot 0.5 mL (minimum 0.4 mL) of serum into a plastic vial and freeze (preferred) or refrigerate up to 5 days.

[View PDF](#)