



## **Test update 870**

**Posted Date** 02/21/2024

**Effective Date** 02/21/2024

**Test Name** [Epstein Barr Virus DNA by PCR, Quantitative, Plasma](#)

**Update Type** [Test Methodology Changed](#)

**CPT Code** 87799

**EFFECTIVE DATE:** February 21, 2024

### **TEST PLATFORM UPDATE**

#### **Epstein Barr Virus DNA by PCR, Quantitative, Plasma**

Order Code: QEBV

Fee Code: LA007 (CPT 87799)

The Michigan Medicine Clinical Microbiology Laboratory has switched from the Abbott M2000 to the Roche Cobas 600 platform for Epstein Barr Virus DNA by PCR, Quantitative, testing effective February 21, 2024.

Please note the range of quantification for the old and new platforms:

#### **Abbott M2000 (old):**

Range of Quantification: 501 - 10,000,000 IU/mL EBV DNA

Range of Quantification: 2.7 - 7.0 Log<sub>10</sub> IU/mL EBV DNA

#### **Roche Cobas 6800 (new):**

Range of Quantification: 35 - 100,000,000 IU/mL EBV DNA

Range of Quantification: 1.54 - 8.0 Log<sub>10</sub> IU/mL EBV DNA

[View PDF](#)