



## **Test Update 901**

**Posted Date** 03/05/2025

**Effective Date** 03/25/2025

**Test Name** [Erythrocytosis Evaluation](#)

**Update Type** [Specimen collection or handling and/or storage changed](#)

**CPT Code** CPT 83020 x2, 83021, 83789

### **SPECIMEN COLLECTION REQUIREMENTS CHANGE**

#### **Erythrocytosis Evaluation**

Order Code: REVE2

Fee Code: AB084, AB085, AA937 (CPT 83020 x2, 83021, 83789)

Reference Laboratory: Mayo REVE2

Effective March 25, 2025, green top (sodium heparin) tubes will no longer be accepted for Mayo Clinic Laboratories Erythrocyte Evaluation assay.

Collection Instructions: Collect specimen in a lavender top (EDTA) (preferred) or yellow top ACD solution B tube. Send 5 mL intact whole blood specimen in original tube. Do not aliquot.

[View PDF](#)