



## **Test Update 683**

**Posted Date** 09/18/2019  
**Effective Date** 09/18/2019  
**Update Type** [New Tests](#)  
**CPT Code** 82017

Order Code: ACARN  
CPT Code: 82017  
Fee Code: DA160

Effective October 1, 2019, MLabs will offer Acylcarnitines, Quantitative, Plasma (ACARN) to replace the Acylcarnitines, Quantitative, Plasma test sent to Mayo Medical Laboratories.

**Specimen Collection:** Collect approximately 3 mL of whole blood in a green top tube (sodium heparin). Centrifuge within 2 hours of collection, aliquot plasma and freeze (preferred) or store at refrigerated temperature (2-8°C) for up to 7 days. Grossly hemolytic specimens are unacceptable. Plasma specimens stored at room temperature greater than 2 hours are unacceptable.

**Methodology:** Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)

**Set Up:** Monday, Thursday

[View PDF](#)