



Test Update 686

Posted Date 10/09/2019

Effective Date 10/15/2019

Update Type [Reference range changed](#)

CPT Code 86225

Order code: ADNA

Effective October 15, 2019, the MLabs will switch the DNA Antibody, Double-Stranded (ADNA) assay from an ELISA format to a chemiluminescent immunoassay methodology. The new assay has fewer false positive results and a better correlation with an accepted reference method Farr assay performed at RDL reference laboratory.

Reference Range

Current Reference Range New Reference Range

< 30.0	< 27 IU/mL Negative
	27-35 IU/mL Borderline Positive
	≥ 36 IU/mL Positive

[View PDF](#)