



Test Update 686

Posted Date 10/09/2019

Effective Date 10/15/2019

Update Type [Reference range changed](#)

CPT Code 86225

Order code: ADNA

Effective October 15, 2019, the MLabs will switch the DNA Antibody, Double-Stranded (ADNA) assay from an ELISA format to a chemiluminescent immunoassay methodology. The new assay has fewer false positive results and a better correlation with an accepted reference method Farr assay performed at RDL reference laboratory.

Reference Range

Current Reference Range New Reference Range

| | |
|--------|---------------------------------|
| < 30.0 | < 27 IU/mL Negative |
| | 27-35 IU/mL Borderline Positive |
| | ≥ 36 IU/mL Positive |

[View PDF](#)