



## [Test Update 682](#)

**Posted Date** 09/13/2019  
**Effective Date** 09/30/2019  
**CPT Code** 87799

Order Code: **QBKV**

Fee Code: LA019

Effective September 30, 2019, MLabs will offer BK Virus DNA, Quantitative, Plasma (QBKV) to replace the ~~V~~**BKP** test sent to Viracor. This test is primarily used for kidney transplant patients to prevent BK virus associated nephropathy and also has a role in bone marrow transplant patients who many develop hemorrhagic cystitis.

**Specimen Collection:** Collect approximately 7 mL of whole blood in a EDTA Lavender Top tube. Centrifuged and aliquot the plasma into a plastic vial. Ship the specimen frozen. This specimen cannot be shared with other testing.

**Methodology:** Polymerase Chain Reaction (PCR)

**Set Up/ Analytic Time:** Monday- Saturday

[View PDF](#)