



## [Test Update 681](#)

**Posted Date** 09/04/2019

**Effective Date** 09/04/2019

**Update Type** [New Tests](#)

**CPT Code** 86038, 86200, if applicable 83516 x2

Order Code: CTDC

CPT Code: 86038, 86200, if applicable 83516 x2

Reference Laboratory: Mayo Medical Laboratory

Effective September 4, 2019, MLabs will offer Connective Tissue Disease Cascade, Serum (CTDC) testing sent to Mayo Medical Laboratories.

Collection Instructions: Collect approximately 3 mL of whole blood in a red top or serum separator tube (SST). Centrifuge and aliquot the serum into a plastic vial. Ship the specimen refrigerated (preferred) or frozen.

Methodology: ELISA

Reference Range:

[View PDF](#)