



Test Update 711

Posted Date 04/01/2020

Effective Date 04/01/2020

Update Type [Specimen collection or handling and/or storage changed](#)

CPT Code 80159

Order Code: CLOZ

CPT Code 80159

Referral Laboratory: Mayo Medical Laboratories

Effective April 1, 2020, MLabs has changed the specimen stability for Clozapine (CLOZ) per Mayo Medical Laboratories.

Current Specimen Stability

Specimen	Temperature	Time
Serum Red	Refrigerated (preferred)	14 Days
	Ambient	14 Days
	Frozen	14 Days

New Specimen Stability

Specimen	Temperature	Time
Serum Red	Refrigerated (preferred)	28 Days
	Ambient	7 Days
	Frozen	28 Days

[View PDF](#)