



## **Test Update 702**

**Posted Date** 02/05/2020  
**Effective Date** 03/03/2020  
**Update Type** [Replacement Test](#)  
**CPT Code** 83864, 82542

### **TEST OBSOLETE**

#### **Mucopolysaccharides Screen, Urine**

**Order Code:** MPSSC  
**CPT Code:** 83864, 82542  
**Reference Laboratory:** Mayo Medical Laboratories

### **REPLACEMENT TEST**

#### **Mucopolysaccharides Quantitative, Random, Urine**

**New Order Code:**  
**CPT Code:** 83864, 82570  
**Reference Laboratory:** Mayo Medical Laboratories

Effective March 3, 2020, Mucopolysaccharides Screen, Urine (MPSSC) will become obsolete and replaced by Mucopolysaccharides Quantitative, Random, Urine (XXXX) per Mayo Medical Laboratories.

Collection Instructions: *Patient Preparation:* Do not administer low-molecular weight heparin prior to collection. Collect approximately 3 mL of random urine preferably in the early morning. Aliquot the specimen in a 5 mL plastic urine tube. Ship the specimen refrigerated.

Methodology: Liquid Chromatography - Tandem Mass Spectrometry (LC-MS/MS)

[View PDF](#)