



## [Test Update 718](#)

**Posted Date** 05/20/2020  
**Effective Date** 06/02/2020  
**Update Type** [New Tests](#)  
**CPT Code** 84163

### **First Trimester Maternal Screen**

Effective June 2, 2020, MLabs will offer First Trimester Maternal Screen (1STT1) testing via Mayo Medical Laboratories.

Collection Instructions: Collect 3 mL of whole blood in a red top vacutainer. Centrifuge and aliquot the serum into a plastic vial within two hours of collection. Ship the specimen refrigerated (preferred) within 7 days, or freeze and ship within 90 days.

Methodology: Immunoenzymatic Assay

Analytic Time: Mon- Fri, 1-3 days

[View PDF](#)