



Test Update 718

Posted Date 05/20/2020
Effective Date 06/02/2020
Update Type [New Tests](#)
CPT Code 84163

First Trimester Maternal Screen

Effective June 2, 2020, MLabs will offer First Trimester Maternal Screen (1STT1) testing via Mayo Medical Laboratories.

Collection Instructions: Collect 3 mL of whole blood in a red top vacutainer. Centrifuge and aliquot the serum into a plastic vial within two hours of collection. Ship the specimen refrigerated (preferred) within 7 days, or freeze and ship within 90 days.

Methodology: Immunoenzymatic Assay

Analytic Time: Mon- Fri, 1-3 days

[View PDF](#)