



Test Update 725

Posted Date 07/29/2020
Effective Date 08/15/2020
Update Type [Replacement Test](#)
CPT Code 80299

TEST CHANGE

TEST OBSOLETE

Warfarin, Serum

Order Code: WRF
CPT Code: 80299
Referral Laboratory: Mayo Laboratories

REPLACEMENT TEST

Warfarin, Plasma

Order Code: FWARP
CPT Code: 80375
Referral Laboratory: Mayo Laboratories

Effective August 15, 2020, MLabs will be replacing the Warfarin, Serum (FSCMS) test and replacing it with Warfarin, Plasma (FWARP) test per Mayo Medical Laboratories.

Collection Instructions: Collect 3mL of whole blood in an EDTA lavender top vacutainer or pink top vacutainer. Centrifuge and aliquot the plasma into a plastic vial. Ship the specimen refrigerated to MLabs.

Specimen Stability

Specimen Type	Temperature	Time
----------------------	--------------------	-------------

Plasma	Refrigerated (preferred)	14 Days
--------	--------------------------	---------

Ambient

14 Days

Frozen

90 Days

Methodology: High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

Test Analysis: Mon- Sun 8 days

Reference Range: 0.2 mcg/mL

The reported value may vary based on dilution or weight of the specimen

[View PDF](#)