



## **Test Update 774**

**Posted Date** 08/25/2021

**Effective Date** 08/16/2021

**Test Name** [Prostate Specific Antigen, Free and Total](#)

**Update Type** [Specimen collection or handling and/or storage changed](#)

**CPT Code** 84153 Total, 84154 Free

### **TEST CHANGE**

#### **Prostate Specific Antigen, Free and Total**

**Order Code:** PSAFT

**CPT Code:** 84153 Total, 84154 Free

**Referral Laboratory:** Mayo Medical Laboratories

Effective August 16, 2021, MLabs has changed the temperature requirements for Prostate Specific Antigen, Free and Total (PSAFT) per Mayo Medical Laboratories.

#### **Current Temperature Requirement**

<b>Specimen</b>	<b>Temperature Time</b>	
Serum	Frozen	90 Days

#### **New Temperature Requirement**

<b>Specimen</b>	<b>Temperature</b>	<b>Time</b>
Serum	Frozen (preferred)	90 Days
	Refrigerated	72 Hours

