



## **Test Update 786**

**Posted Date** 12/02/2021

**Effective Date** 12/01/2021

**Test Name** [Epstein Barr Viral Capsid Antigen, IgM](#)

**Update Type** [Test Down or Delayed](#)

**CPT Code** 86665

### **TEST DOWN**

#### **Epstein Barr Viral Capsid Antigen, IgM**

**Order Code:** EBVCM

**CPT Code:** 86665

Effective December 1, 2021, MLabs EBV Viral Capsid Ag IgM (EBVCM) assay is down effective 12/1/21 due to reagent backorder. Requests for Epstein Barr Viral Capsid Antigen, IgM (EBVCM) will be sent to Warde until further notice.

Collection Instructions: Collect approximately 2 mL of whole blood in a serum separator tube (SST) or red top vacutainer, centrifuge and aliquot the serum into a plastic vial. Ship specimen refrigerated.

[View PDF](#)