



## **Test Update 788**

**Posted Date** 12/15/2021  
**Effective Date** 12/14/2021  
**Test Name** [Pancreatic Polypeptide](#)  
**Update Type** [Replacement Test](#)  
**CPT Code** 82397

### **TEST DELAY**

#### **Pancreatic Polypeptide**

**Order Code:** HPP  
**CPT Code:** 83519  
**Referral Laboratory:** Mayo Medical Laboratories

Effective December 14, 2021 Pancreatic Polypeptide (HPP) test will be delayed due to performance issues. Specimens will be stabilized and tested when the issue is resolved per Mayo Medical Laboratories.

### **RECOMMENDED ALTERNATIVE**

#### **Pancreatic Polypeptide, ECL**

**Order Code:** MMLR  
**CPT Code:** 82397  
**Referral Laboratory:** Mayo Medical Laboratories

Collection Instructions: Collect approximately 6 mL of blood in a lavender top vacutainer. Centrifuge and aliquot the plasma in a plastic tube. Ship the specimen frozen, ambient or refrigerated.

Methodology: Electrochemiluminescence (ECLIA)

Analytic Time: 4-8 days, test performed on Wednesday

**Reference Range:** 200-1000 ng/L • Ann Arbor, MI 48109-2800 • 800.862.7284 • 734.936.0755 Fax • [mlabs.umich.edu](http://mlabs.umich.edu)

<b>Age</b>	<b>Reference Range</b>
------------	------------------------

<4 years	Not established
----------	-----------------

4-17 years	92-752 pg/mL
------------	--------------

≥18 years	56-480 pg/mL
-----------	--------------

[View PDF](#)