



Michigan Medicine  
Laboratories (MLabs)

N-LNC Specimen Processing  
2800PlymouthRd,Bldg35  
AnnArbor,MI48109-2800

734.936.2598 • 800.862.7284 • mlabs.umich.edu

FAX: 734.936.0755

# PHYSICIAN OFFICE SUPPLY REQUISITION

Date Requested: \_\_\_\_\_

Client Location Code: \_\_\_\_\_

Client Name: \_\_\_\_\_

Contact Name: \_\_\_\_\_

Phone Number: \_\_\_\_\_

MLabs provides selected supplies only for the collection and transport of specimens to MLabs. The number of items ordered must reasonably correlate with the number of specimens sent to MLabs and will be monitored and adjusted as needed.

Stock #	BLOOD COLLECTION	Circle Quantity	Order	Stock #	MICROBIOLOGY/VIROLOGY	Circle Quantity	Order
PB 62	Gold SST 3.5 mL	Bx 100 or Each		PM 5	Aerobic Culture Swab (ESwab)	Bag 50 or Each	
PB 16	Gold SST 5 mL	Bx 100 or Each		PM 6	Anaerobic Culture Transport Tube	Each	
PB 21	Red 4 mL (no additive)	Bx 100 or Each		PM 15	CT/GC/TV/MG – Aptima Unisex Swab Kit Endcx/Urethra/Oral/Rectal	Each	
PB 15	Red 6 mL (no additive)	Bx 100 or Each		PM 20	CT/GC/TV/MG – Aptima Multitest (Vag/Rectal) Kit	Each	
PB 56	Lavender EDTA 3 mL (use for CBC)	Bx 100 or Each		PM 18	CT/GC/TV/MG – Aptima Urine Tube/Pipette	Each	
PB 42	Lavender EDTA 6 mL	Bx 100 or Each		PM 35	M4-RT Transport Kits (w/swab)	Each	
PB 12	Lt. Blue w/sodium citrate 2.7 mL	Bx 100 or Each		PM 19	Mini-tip flocced Swabs (Nasopharyngeal)	Each	
PB 36	Tan w/EDTA 3 mL (Lead Testing)	Each		PM 21	Blood Culture Aerobic Bottle (Blue Top)	Each	
PB 53	Green (sodium heparin) 4 mL	Each		PM 23	Blood Culture Anaerobic Bottle (Purple Top)	Each	
PB 52	Green (sodium heparin) 6 mL	Each					
PB 46	Royal (Navy) Blue, no additive 7 mL	Each					
PB 75	Royal (Navy) Blue, w/EDTA 7 mL	Each					
PB 65	Yellow ACD Solution A 8.5 mL	Each		<b>FORMS</b>			
PB 68	Pink (EDTA-Blood Bank) 6 mL	Each		MLO	Physician Office Requisition	Each	
PB 30	Straight Needles 21g, 1 1/4"	Box of 48		MLO	Dermatopathology Requisition	Each	
PB 67	Straight Needles 22g, 1"	Box of 48		MLO	Surgical/Cytopathology Requisition	Each	
PB 29A	Butterfly Needles with safety lock 21g	Box of 50		MLO	Allergen Requisition	Each	
PB 29B	Butterfly Needles with safety lock 23g	Box of 50		MLO	MMGL BRCA Requisition	Each	
PB 28A	Tourniquets (Latex Free)	Pack of 50		MLO	QUAD/AFP Requisition	Each	
PB 24	Vacutainer Holder	Bag of 250		MLO	Physician Office Supply Requisition	Each	
PX 10	Glucola: 50g	Box 6 or Each		MLO	Drawing Station Maps	Pack of 50	
PX 190	Glucola: 75g	Box 6 or Each		MLO	Lead Forms	Each	
PX 11	Glucola: 100g	Box 6 or Each		MLO	HIV Consent Forms	Pack of 25	
				MLO	Stool Specimen Collection Instructions	Pack of 50	
<b>URINE / STOOL COLLECTION</b>				<b>MISCELLANEOUS</b>			
PX 14	Urine Cup, Sterile	Bag of 75		PX 70	Specimen Bags, Refrigerated	Pack of 100	
PB 79	Urine Culture Vacutainer Transport (Gray Top)	Each		PX 71	Specimen Bags, Room Temperature	Pack of 100	
PX 5	24 Hour Urine Jug, Plain	Each		PX 69	Specimen Bags, Red STAT	Pack of 10	
PX 24	24 Hour Urine Jug w/Sodium Carbonate*	Each		PX 59	Large Ziploc Biohazard Bag (11" x 18")	Each	
PX 240	24 Hour Urine Jug w/Glacial Acetic Acid*	Each		PF 20	Specimen Labels (2" x 1 1/4")	Roll of 100	
PX 530	24 Hour Fecal Fat Canister	Each		PX 505	Sharps Container, Small	Each	
PX 34	Fecal Occult Bld FOBT Bottle (on-site collect)	Each or Bag 50		PF 45	MLabs Connect Barcode Labels	Roll	
PX 55	Fecal Occult Bld FIT kit w/mailler (home collect)	Each or Box 20					
PM 4	Ova & Parasite Kit	Each		<b>SURGICAL PATHOLOGY</b>			
PM 30	Pinworm Collection Tube*	Each		PX 210	Formalin 10 mL (20 mL vial)	Box 24 or Each	
PM 12	Stool Culture Kit	Each		PX 41	Formalin 60 mL vial	Box 24 or Each	
				IM 380	Zeus Tissue Fixative 10 mL vial*	Each	
<b>CYTOLOGY</b>				<b>OTHER</b>			
PK 36	Thin Prep/HPV Kits (vial, brush & spatula included)*	Box of 25 Kits					
PX 63	Cytology Brush (Not for ThinPrep Use)	Box of 100					
CY140	CytoLyt Solution 30 mL (Fine Needle Aspirate)	Each					

\* Item supplied by Pathology Stores/Material Services

White = MLabs    Yellow = Client

Revised: 02-05-20